


Product Datasheet

RRM2B antibody (orb536297)

Description	Rabbit polyclonal antibody to RRM2B
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human RRM2B(NP_056528.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	4kDa/7kDa/15kDa/34kDa/40kDa/48kDa
Uniprot ID	Q7LG56
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A549 cell...