



Product Datasheet

Fibrinogen antibody (orb536312)


Description	Rabbit Polyclonal antibody to Fibrinogen
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human Fibrinogen alpha chain (FGA)(NP_000499.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	95KDa/60-70KDa
Uniprot ID	P02671
Dilution Range	IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U-2 OS cells...



Immunofluorescence analysis of NIH/3T3 cells...



Immunofluorescence analysis of C6 cells...