



Product Datasheet

SLC25A27 antibody (orb536319)


Description	Rabbit polyclonal antibody to SLC25A27
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic peptide of human SLC25A27(NP_004268.3).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36kDa
Uniprot ID	O95847
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH/3T3 cells...



Immunofluorescence analysis of HeLa cells...



Immunofluorescence analysis of C6 cells...