


Product Datasheet

AP2M1 antibody (orb536331)

Description	Rabbit polyclonal antibody to AP2M1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human AP2M1(NP_004059.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50KDa
Uniprot ID	Q96CW1
Dilution Range	IF: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...