


Product Datasheet

HRNR antibody (orb536367)

Description	Rabbit polyclonal antibody to HRNR
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human HRNR(NP_001009931.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	282kDa
Uniprot ID	Q86YZ3
Dilution Range	IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...