


Product Datasheet

SLC24A5 antibody (orb536385)

Description	Rabbit polyclonal antibody to SLC24A5
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human SLC24A5(NP_995322.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48kDa/54kDa
Uniprot ID	Q71RS6
Dilution Range	IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...