


Product Datasheet

ST6GAL1 antibody (orb536391)

Description	Rabbit polyclonal antibody to ST6GAL1
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human ST6GAL1(NP_775323.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20kDa/46kDa
Uniprot ID	P15907
Dilution Range	IF: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...