

Product Datasheet

EGLN1 antibody (orb536405)

Description	Rabbit polyclonal antibody to EGLN1
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human PHD2/EGLN1(NP_071334.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36kDa/43kDa/46kDa
Uniprot ID	Q9GZT9
Dilution Range	IF: 1:50-1:200
Expiration Date	12 months from date of receipt.