

## Product Datasheet

### ABL1 (Phospho-Y245) antibody (orb536408)

|                            |                                                                            |
|----------------------------|----------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to ABL1 (Phospho-Y245)                          |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                          |
| <b>Conjugation</b>         | Unconjugated                                                               |
| <b>Tested Applications</b> | IHC, WB                                                                    |
| <b>Immunogen</b>           | A synthetic phosphorylated peptide around Y245 of human ABL1(NP_005148.2). |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.                           |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                                |
| <b>Note</b>                | For research use only                                                      |
| <b>Isotype</b>             | IgG                                                                        |
| <b>MW</b>                  | 124kDa                                                                     |
| <b>Uniprot ID</b>          | <b>P00519</b>                                                              |
| <b>Dilution Range</b>      | WB: 1:500-1:2000, IF: 1:50-1:200                                           |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                            |