


Product Datasheet

Dnmt1 (Phospho-pY399) antibody (orb536410)

Description	Rabbit polyclonal antibody to Dnmt1 (Phospho-pY399)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y399 of human Phospho-DNMT1-pY399 (NP_001370.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P26358
Entrez	1786
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MCF-7 cel...