


Product Datasheet

Neuropeptide Y antibody (orb536436)

Description	Rabbit Polyclonal antibody to Neuropeptide Y
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human NPY(NP_000896.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	10kDa
Uniprot ID	P01303
Dilution Range	IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U-251MG C...