


Product Datasheet

CYP11B1 antibody (orb536456)

Description	Rabbit polyclonal antibody to CYP11B1
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC
Immunogen	Recombinant fusion protein of human CYP11B1(NP_000488.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	49kDa/57kDa
Uniprot ID	P15538
Dilution Range	IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...