

Product Datasheet

Cd207 antibody (Biotin) (orb53646)

Description	Cd207 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Mouse C-type lectin domain family 4 member K protein (63-331AA)
Target	Cd207
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8VBX4
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.