


Product Datasheet

ENTPD5 antibody (orb536470)

Description	Rabbit polyclonal antibody to ENTPD5
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human ENTPD5(NP_001240.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	47kDa
Uniprot ID	075356
Dilution Range	IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...