


Product Datasheet

MAP4K5 antibody (orb536492)

Description	Rabbit polyclonal antibody to MAP4K5
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human MAP4K5(NP_006566.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	95kDa
Uniprot ID	Q9Y4K4
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...