


Product Datasheet

SH2D1B antibody (orb536525)

Description	Rabbit polyclonal antibody to SH2D1B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Recombinant fusion protein of human SH2D1B (NP_444512.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	O14796
Entrez	117157
Dilution Range	IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence
analysis of U2OS
cell...