



Product Datasheet

YBX1 antibody (orb536552)


Description	Rabbit polyclonal antibody to YBX1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic peptide of human YB-1/YB-1/YBX1(NP_004550.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	49kDa
Uniprot ID	P67809
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of PC-12 cell...



Immunofluorescence analysis of NIH/3T3 C...



Immunofluorescence analysis of HeLa cell...