



Product Datasheet

ARSB antibody (orb536565)


Description	Rabbit polyclonal antibody to ARSB
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	Recombinant fusion protein of human ARSB (NP_000037.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P15848
Entrez	411
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse brain tissue.



Immunohistochemical staining of mouse liver tissue.



Immunohistochemical staining of human tissue.