


Product Datasheet

PRKAG3 antibody (orb536580)

Description	Rabbit polyclonal antibody to PRKAG3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Recombinant fusion protein of human PRKAG3 (NP_059127.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q9UGI9
Entrez	53632
Dilution Range	IHC-P: 1:50-1:200, IF: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human kidney tissue.