



Product Datasheet

ATP5I antibody (orb536591)


Description	Rabbit polyclonal antibody to ATP5I
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human ATP5I(NP_009031.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	7kDa
Uniprot ID	P56385
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ut...



Immunohistochemical staining of human st...



Immunohistochemical staining of human br...