



Product Datasheet

ITGB3BP antibody (orb536595)


Description	Rabbit polyclonal antibody to ITGB3BP
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human ITGB3BP(NP_055103.3).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	6kDa/12kDa/19kDa/20kDa/24kDa
Uniprot ID	Q13352
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse brain tissue.



Immunohistochemical staining of human gastric tissue.



Immunohistochemical staining of human appendix tissue.