



Product Datasheet

SLC5A6 antibody (orb536600)


Description	Rabbit polyclonal antibody to SLC5A6
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	A synthetic peptide of human SLC5A6(NP_066918.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	68kDa
Uniprot ID	Q9Y289
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse heart tissue.



Immunohistochemical staining of rat heart tissue.



Immunohistochemical staining of human liver tissue.