



Product Datasheet

GAB1 antibody (orb536606)

Description	Rabbit polyclonal antibody to GAB1
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human GAB1(NP_997006.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	76kDa/80kDa
Uniprot ID	Q13480
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human placenta.



Immunohistochemical staining of human lung tissue.