

Product Datasheet

DME antibody (FITC) (orb53661)

Description	DME antibody (FITC)
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	FITC
Immunogen	Recombinant Arabidopsis thaliana Transcriptional activator DEMETER protein (1-320AA)
Target	DME
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8LK56
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.