



Product Datasheet

TAX1BP1 antibody (orb536610)

Description	Rabbit polyclonal antibody to TAX1BP1
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human TRAF6BP/TRAF6BP/TAX1BP1(NP_006015.4).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65kDa/68kDa/86kDa/90kDa
Uniprot ID	Q86VP1
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human brain tissue.



Immunohistochemical staining of human lung tissue.