



Product Datasheet

MTA2 antibody (orb536639)


Description	Rabbit polyclonal antibody to MTA2
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human MTA2(NP_004730.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	55kDa/75kDa
Uniprot ID	O94776
Dilution Range	WB: 1:200-1:1000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human gastric tissue.



Immunohistochemical staining of human lung tissue.



Immunohistochemical staining of human tissue.