



Product Datasheet

THBS1 antibody (orb536665)


Description	Rabbit polyclonal antibody to THBS1
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	A synthetic peptide of human THBS1(NP_003237.2).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	120kDa/129kDa
Uniprot ID	P07996
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse spleen



Immunohistochemical staining of mouse brain



Immunohistochemical staining of mouse spleen