



Product Datasheet

DDX21 antibody (orb536668)


Description	Rabbit polyclonal antibody to DDX21
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC
Immunogen	Recombinant fusion protein of human DDX21(NP_004719.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	79kDa/87kDa
Uniprot ID	Q9NR30
Dilution Range	IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U-2 OS cells.



Immunohistochemical staining of mouse brain tissue.



Immunohistochemical staining of mouse testis tissue.