



Product Datasheet

CYP17A1 antibody (orb536672)

Description	Rabbit polyclonal antibody to CYP17A1
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CYP17A1(NP_000093.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	57kDa
Uniprot ID	P05093
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat brain...



Western blot analysis of extracts of rat...