



Product Datasheet

GRIN1 antibody (orb536675)

Description	Rabbit polyclonal antibody to GRIN1
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	Recombinant fusion protein of human NMDAR1(NP_001172019.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	99kDa/101kDa/103kDa/105kDa/106kDa/107kDa
Uniprot ID	Q05586
Dilution Range	IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...



Immunohistochemical staining of rat brain...