



Product Datasheet

WNT3A antibody (orb536699)


Description	Rabbit polyclonal antibody to WNT3A
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human WNT3A(NP_149122.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	39kDa/42kDa
Uniprot ID	P56704
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human gastric tissue.



Immunohistochemical staining of human lung tissue.



Immunohistochemical staining of rat brain tissue.