



Product Datasheet

LRRK2 antibody (orb536722)


Description	Rabbit polyclonal antibody to LRRK2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human LRRK2(NP_940980.3).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	290KDa
Uniprot ID	Q5S007
Dilution Range	IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse spinal cord tissue.



Immunohistochemical staining of mouse lung tissue.



Immunohistochemical staining of rat heart tissue.