



Product Datasheet

ZNF177 antibody (orb536747)


Description	Rabbit polyclonal antibody to ZNF177
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC
Immunogen	Recombinant fusion protein of human ZNF177(NP_003442.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	13kDa/36kDa/54kDa
Uniprot ID	Q13360
Dilution Range	IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of C6 cells ...



Immunohistochemical staining of mouse pa...



Immunohistochemical staining of mouse ki...