

Product Datasheet

ATRIP antibody (orb536764)

Description	Rabbit polyclonal antibody to ATRIP
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human ATRIP (NP_569055.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	100kDa
Uniprot ID	Q8WXE1
Entrez	84126
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200, IP: 1:20-1:50
Expiration Date	12 months from date of receipt.