

Product Datasheet

Lck (Phospho-Y505) antibody (orb536765)

Description	Rabbit polyclonal antibody to Lck (Phospho-Y505)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic phosphorylated peptide around Y505 of human Lck(NP_001036236.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	56KDa
Uniprot ID	P06239
Dilution Range	WB: 1:500-1:2000, IP: 1:50-1:100
Expiration Date	12 months from date of receipt.