



Product Datasheet

PTPN11 (Phospho-Y542) antibody (orb536772)

Description	Rabbit polyclonal antibody to PTPN11 (Phospho-Y542)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic phosphorylated peptide around Y542 of human PTPN11(NP_002825.3).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	72kDa
Uniprot ID	Q06124
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human brain tissue.



Western blot analysis of extracts from 293T cells.