


Product Datasheet

RELA (Phospho-S536) antibody (orb536781)

Description	Rabbit polyclonal antibody to RELA (Phospho-S536)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S536 of human RelA(NP_068810.3).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	65KDa
Uniprot ID	Q04206
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100, IF: 1:100-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from H...