


Product Datasheet

EIF4EBP1 (Phospho-T36) antibody (orb536787)

Description	Rabbit polyclonal antibody to EIF4EBP1 (Phospho-T36)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding T36 of human Phospho-eIF4EBP1-T36
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	20kDa
Uniprot ID	Q13541
Entrez	1978
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from M...