


Product Datasheet

DPP6 antibody (orb536833)

Description	Rabbit polyclonal antibody to DPP6
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human DPP6(NP_001277182.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	105kDa
Uniprot ID	P42658
Dilution Range	WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of 293...