

Product Datasheet

PIGK antibody (HRP) (orb53684)

Description	PIGK antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human GPI-anchor transamidase protein (91-293AA)
Target	PIGK
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q92643
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.