


Product Datasheet

IGF1R antibody (orb536846)

Description	Rabbit polyclonal antibody to IGF1R
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human IGF1R(NP_000866.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	100kDa
Uniprot ID	P08069
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:100, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...