



Product Datasheet

KAT6A antibody (orb536849)


Description	Rabbit polyclonal antibody to KAT6A
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human KAT6A(NP_001092883.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	250kDa
Uniprot ID	Q92794
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of 293...



Immunofluorescence analysis of NIH/3T3 C...



Immunofluorescence analysis of U2OS cell...