



Product Datasheet

ATM (Phospho-S1981) antibody (orb536873)


Description	Rabbit polyclonal antibody to ATM (Phospho-S1981)
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S1981 of human ATM(NP_000042.3).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	370KDa
Uniprot ID	Q13315
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...



Immunofluorescence analysis of U2OS cell...



Immunofluorescence analysis of MCF-7 cel...