


Product Datasheet

VASP (Phospho-S239) antibody (orb536884)

Description	Rabbit polyclonal antibody to VASP (Phospho-S239)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic phosphorylated peptide around S239 of human VASP(NP_003361.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50KDa
Uniprot ID	P50552
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of 293...