



Product Datasheet

SLC4A4 antibody (orb536900)


Description	Rabbit polyclonal antibody to SLC4A4
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic peptide of human SLC4A4 / NBC(NP_001091954.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	135KDa
Uniprot ID	Q9Y6R1
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse kidney tissue.



Immunohistochemical staining of rat kidney tissue.



Immunohistochemical staining of human prostate tissue.