


Product Datasheet

EGFR (Phospho-Y1045) antibody (orb536938)

Description	Rabbit polyclonal antibody to EGFR (Phospho-Y1045)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y1045 of human EGFR(NP_005219.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	180kDa
Uniprot ID	P00533
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of A-4...