


Product Datasheet

STK11 (Phospho-T189) antibody (orb536941)

Description	Rabbit polyclonal antibody to STK11 (Phospho-T189)
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around T189 of human STK11(NP_000446.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	54kDa
Uniprot ID	Q15831
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of A431...