



## Product Datasheet

### PAI-1 antibody (orb536973)

|                            |                                                                    |
|----------------------------|--------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to SERPINE1                             |
| <b>Species/Host</b>        | Rabbit                                                             |
| <b>Reactivity</b>          | Human, Mouse                                                       |
| <b>Conjugation</b>         | Unconjugated                                                       |
| <b>Tested Applications</b> | IF, WB                                                             |
| <b>Immunogen</b>           | Recombinant fusion protein of human PAI-1/Serpin E1 (NP_000593.1). |
| <b>Preservatives</b>       | Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.         |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                        |
| <b>Note</b>                | For research use only                                              |
| <b>Isotype</b>             | IgG                                                                |
| <b>Clonality</b>           | Polyclonal                                                         |
| <b>Purity</b>              | Affinity purification                                              |
| <b>MW</b>                  | 45kDa                                                              |
| <b>Uniprot ID</b>          | <a href="#">P05121</a>                                             |
| <b>Entrez</b>              | <a href="#">5054</a>                                               |
| <b>Dilution Range</b>      | WB: 1:500-1:2000, IF: 1:50-1:200                                   |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                    |



Immunofluorescence analysis of human liv...



Western blot analysis of extracts of A-5...