


Product Datasheet

SNX17 antibody (orb536976)

Description	Rabbit polyclonal antibody to SNX17
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SNX17(NP_055563.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45kDa
Uniprot ID	Q15036
Dilution Range	WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of A-549 cells.