


Product Datasheet

TBX21 antibody (orb536988)

Description	Rabbit polyclonal antibody to TBX21
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human TBX21(NP_037483.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	58KDa
Uniprot ID	Q9UL17
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of B-
C...